
# eFlow<sup>®</sup> rapid nebuliser system

CE 0123

eFlow<sup>®</sup>  
TECHNOLOGY

**en** Instructions for Use

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**sk** Návod na použitie

**sl** Navodila za uporabo

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**Patented.** See [www.pari.com/ip](http://www.pari.com/ip)

©2023 PARI Pharma GmbH. All rights reserved.

The PARI Pharma logo is a registered trademark of PARI GmbH. eFlow<sup>®</sup>, Tolero<sup>®</sup> and Vantobra<sup>®</sup> are registered trademarks of PARI Pharma GmbH. TouchSpray made under licence from The Technology Partnership plc.

<b>en</b>	<b>INSTRUCTIONS FOR USE</b> .....	<b>3</b>
<b>bg</b>	<b>ИНСТРУКЦИИ ЗА УПОТРЕБА</b> .....	<b>43</b>
<b>cs</b>	<b>NÁVOD K POUŽITÍ</b> .....	<b>85</b>
<b>el</b>	<b>ΟΔΗΓΙΕΣ ΧΡΗΣΕΩΣ</b> .....	<b>125</b>
<b>et</b>	<b>KASUTUSJUHEND</b> .....	<b>167</b>
<b>hr</b>	<b>UPUTE ZA UPORABU</b> .....	<b>207</b>
<b>hu</b>	<b>HASZNÁLATI UTASÍTÁS</b> .....	<b>247</b>
<b>is</b>	<b>NOTKUNARLEIÐBEININGAR</b> .....	<b>287</b>
<b>pl</b>	<b>INSTRUKCJA OBSŁUGI</b> .....	<b>327</b>
<b>sk</b>	<b>NÁVOD NA POUŽITIE</b> .....	<b>367</b>
<b>sl</b>	<b>NAVODILA ZA UPORABO</b> .....	<b>407</b>
<b>tr</b>	<b>KULLANIM TALİMATLARI</b> .....	<b>447</b>

# en INSTRUCTIONS FOR USE

## For eFlow<sup>®</sup>rapid nebuliser system 178G1102

### Table of contents

<b>A Introduction</b> .....	<b>4</b>
<b>B Safety Precautions</b> .....	<b>6</b>
<b>C Parts of Your eFlow<sup>®</sup>rapid nebuliser system</b> .....	<b>8</b>
<b>D Charging Your Controller</b> .....	<b>10</b>
<b>E Cleaning Your Nebuliser Handset</b> .....	<b>12</b>
<b>F Disinfecting Your Nebuliser Handset</b> .....	<b>15</b>
<b>G Cleaning and Disinfecting Your Controller, Connection Cord, Power Supply</b> .....	<b>18</b>
<b>H Assembling Your Nebuliser System</b> .....	<b>19</b>
<b>I Taking Your Treatment</b> .....	<b>22</b>
<b>J Controller Settings</b> .....	<b>26</b>
<b>K Troubleshooting</b> .....	<b>29</b>
<b>L Terms and Conditions of Warranty</b> .....	<b>33</b>
<b>M Spare Parts and Accessories</b> .....	<b>34</b>
<b>N Technical Data</b> .....	<b>35</b>
<b>O Symbols</b> .....	<b>38</b>
<b>Quick Instructions</b> .....	<b>40</b>

## A Introduction


The eFlow<sup>®</sup>*rapid* nebuliser system is a quiet, lightweight, battery-operated medical device which converts your drug into an aerosol for inhalation.

The drug is placed into the medication reservoir of your nebuliser handset, which feeds to a membrane containing tiny holes. When the device is turned on, the membrane vibrates forcing the drug through the tiny holes to form a fine aerosol mist, which you breathe into your lungs.

The eBase<sup>®</sup> Controller is used to operate eFlow<sup>®</sup> Technology nebuliser handsets, e.g. Tolero<sup>®</sup> for Vantobra<sup>®</sup> and eFlow<sup>®</sup>*rapid*.

This IFU depicts the eFlow<sup>®</sup>*rapid* nebuliser handset and, if not explicitly described, the instructions must also be followed when using a drug-specific nebuliser handset.

The drug-specific nebuliser handset is made available with a specially approved drug and/or via a specific distribution channel.

This **"Instructions For Use" (IFU)** contains information and safety precautions for your nebuliser system.

Take note of any limitations in the package leaflet accompanying the medication.

Before using your eFlow<sup>®</sup>*rapid* nebuliser system, you must read and understand the entire **"Instructions For Use"**, and save it for future reference.

**If for any reason you do not understand any part of these instructions, contact your service partner before proceeding with your treatment.**

Take special note of all safety precautions marked **Warning or Caution**.

Within the **European Union (EU)**, report serious incidents to the manufacturer and the competent authority.

If your health condition for any reason is not improved, or even worsens, seek professional medical advice.

### Indications

The eFlow<sup>®</sup>*rapid* nebuliser system is indicated to administer solutions or suspensions for a nebuliser handset which are approved for the treatment of diseases of the lower airways.

### Contraindications

No contraindication is known to the manufacturer of this medical device.

### Intended Purpose

The eFlow<sup>®</sup>*rapid* nebuliser system is intended to administer solutions or suspensions for a nebuliser handset as an aerosol for inhalation.

## Intended use

The eFlow<sup>®</sup>*rapid* nebuliser system must only be used in accordance with its intended purpose.

## Intended Patient Group

The eFlow<sup>®</sup>*rapid* nebuliser system is suitable for all age groups but is limited to the patient group defined for the drug. The nebuliser handset is only suitable for persons who are conscious and breathe spontaneously.

## Accessories

For babies and infants who cannot use the mouthpiece to inhale, PARI offers appropriate masks. A mouthpiece with filter reduces the release of the exhaled aerosol into the ambient air. To clean the aerosol head by backwashing the membrane, an *easycare* cleaning aid can be used. Accessories are not available in all countries/sales regions. Contact your service partner if you are unsure whether you should use accessories. When using accessories, follow the instructions that came with the accessories.

## B Safety Precautions

**Please read the following before using this device.**

**Warning** indicates hazards which, without precautionary measures, can result in serious injury or death.

**Caution** indicates hazards which, without precautionary measures, can result in minor or moderate injury or impair treatment; or cause damage to the device.

### WARNING

To reduce the risk of serious injury:

1. **Do not** use a damaged power supply or handle a power supply with wet hands.
2. **Do not** use extension cords where they can fall into a bathtub, sink, or other liquid.
3. Also, **do not** use the nebuliser system where it can fall into water, e.g. bathtub.
4. **Do not** reach for the device if it has fallen into water or other liquid. Unplug immediately from the wall outlet. Retrieve the device **only** after it has been unplugged.
5. Improper use may cause damage to the product and lead to injuries and illness. Follow these Instructions for Use and all safety instructions contained herein and then save them for future reference.
6. To ensure the correct dose, **never** use another drug for a drug-specific nebuliser handset than the one approved, and **do not** use with an eFlow<sup>®</sup>*rapid* nebuliser handset.
7. For reasons of hygiene and to avoid cross-infections, the eFlow<sup>®</sup>*rapid* nebuliser handset and the connection cord must be used by one patient only. For patients at greater risk of infection, e.g., with cystic fibrosis (CF) or a weakened immune system, or for those infected with MRSA it is recommended that the eBase<sup>®</sup> Controller and the power adapter are also used by one patient only.
8. To avoid infections, clean the nebuliser handset including the aerosol head before they are used for the first time.
9. To keep the membrane from clogging and to prevent infections, you must clean the nebuliser handset including the aerosol head immediately after each treatment.
10. To prevent possible infections, you must also disinfect the nebuliser handset and aerosol head after each treatment day.
11. Always unplug the device right after using and before cleaning.
12. To avoid the chance of fire, **do not** use where oxygen or combustible gases are being used such as an oxygen tent.
13. Use only original PARI parts and accessories. Otherwise this can lead to malfunction including interference of nearby devices.
14. The device may be operated in a car, train or airplane. To minimize possible interference issues, please note the following:
  - a. It must only be used in passenger areas of trains and airplanes.
  - b. **Do not** use near the control systems of airplanes or trains.
  - c. It must be operated using battery power while in a car.
15. **Do not** use the device while driving.
16. **Do not** place the controller next to or on top of other medical devices during operation unless both devices are monitored constantly to make sure both are operating properly.
17. **Do not** use in areas exposed to elevated magnetic or electrical radiation such as an MRI scanner or high frequency surgical equipment in hospital.

18. Close supervision is necessary when this product is used by or near children or the physically or mentally impaired persons.
19. Small or broken parts can be swallowed or block the airways and the connection cord poses a strangulation hazard.
20. **Never** operate the nebuliser system if parts or cables appear damaged, deformed or heavily discolored.

 **CAUTION**


To minimize risk of damaging the device or impairing treatment:

1. The controller electronics can be damaged if liquid gets into the controller: **Do not** use the device if liquid gets inside the controller.
2. **Do not** disassemble the controller or modify the equipment at any time. There are no serviceable parts.
3. **Do not** allow pets near the nebuliser system.
4. Portable wireless communication devices (such as cell phones or external antennas) may reduce the performance of the nebuliser system. Keep a minimum distance of 30 cm or 12 inches of these devices and the nebuliser system.
5. Anti-theft systems and Radio Frequency Identification (RFID) readers may reduce performance of the nebuliser system: **Do not** use near entrances of stores, libraries, and hospitals.


## C Parts of Your eFlow<sup>®</sup>rapid nebuliser system

Make sure you have the parts shown below. Several parts are delivered in a carrying case which can also be used to transport the system. The nebuliser system contains an *easycare* cleaning aid.


### Nebuliser Handset (2x)


Aerosol Head


Nebuliser Handset Body


Medication Cap


Blue Inspiratory Valve


Mouthpiece with blue Expiratory Valve

eBase<sup>®</sup> Controller


Connection Cord

Power Supply  
(EU shown)


## General Information About the Nebuliser Handset and Controller

The nebuliser handset body is pre-assembled when delivered and contains the blue inspiratory valve.


You can open the nebuliser handset body by gently pulling up on the tab (c).

The medication reservoir (a) and the aerosol chamber (b) are attached using a "hinge" and must not be detached from one another.

### ⚠ CAUTION

**Do not** touch the silver center part of the aerosol head. Touching it may damage the membrane that generates the aerosol.

The ON/OFF button (a) in the top right corner of the device turns controller on and off.


## D Charging Your Controller

The controller runs with an integrated battery pack and can be recharged at any time.

**Using the power supply:** The power supply will automatically adjust to the incoming voltage and will charge the controller's integrated battery pack. It can be used worldwide, but may require changing the country-specific adapter.

**Step D1** To change the country-specific adapter (a), press the clip (b) and slide the country-specific adapter upwards. Slide on the new country specific adapter until the clip connects.


Fig. D1

**Step D2** Plug the power supply into the controller.


The power cord socket is located on the top of the controller. Push the end of the power supply cord into the power cord socket (a). Do not attempt to insert the power cord into the bottom of the controller.

Fig. D2

**Step D3 Plug the power supply into the wall outlet.**

The power supply will charge the batteries in the controller. As soon as the charging process is completed, you will see a fully charged battery.

*Fig. D3*

**Note:** To preserve battery life, it is recommended to:

- charge battery only when nearly empty.
- unplug the power supply once the battery is fully charged.
- avoid higher temperatures than room temperature during charging.

## E Cleaning Your Nebuliser Handset

### ⚠ WARNING

Your new nebuliser handset is not ready for use. You must clean your nebuliser handset before you use it for the first time and after every use. You must also disinfect your nebuliser handset at the end of the treatment day.

### ⚠ CAUTION

**Do not** touch the silver center part of the aerosol head. Touching it may damage the membrane that generates the aerosol.

**Step E1**     **Disassemble the nebuliser handset and rinse the aerosol head** for about 5 seconds on each side under warm running tap water (quality of drinking water; about 37 °C or 98 °F, warm to the touch).

**Note:** Thorough rinsing of both sides of the aerosol head helps to prevent clogging. Clogging may increase nebulisation time. It is important to clean the aerosol head immediately after each use.

**Step E2**     **Rinse the opened nebuliser handset body (a), the mouthpiece (b), the blue inspiratory valve (c) and the medication cap (d)** under warm running tap water for about 5 seconds to remove medication residue.

Fig. E1



Fig. E2


- Step E3** Add a few drops of **clear liquid dish soap** to warm (about 37 °C or 98 °F) tap water in a clean tub or bowl and **soak all parts of the nebuliser handset for 5 minutes.**

**⚠ CAUTION**

The membrane of the aerosol head can be damaged very easily.


- **Do not** attempt to clean the aerosol head with a brush.
- **Do not** attempt to clean or disinfect the aerosol head in a home dishwasher.
- **Do not** attempt to disinfect the aerosol head in a microwave.

Fig. E3


- Step E4** **Rinse all parts thoroughly** under warm running tap water and then inspect them to make sure they are clean and not visibly damaged. If parts do not appear to be clean, soak them for another 5 minutes as described in Step E3. **Do not** use dirty or damaged parts.

Fig. E4

- Step E5** **Shake** the parts to remove excess water. **Air-dry parts** on a clean, lint-free towel. After the parts are completely dry, store them in a clean bag or cover them with a clean, lint-free towel. Reassemble just before treatment.

Fig. E5

**Note:** If necessary, clean the membrane of the aerosol head with the *easycare* cleaning aid (not more than twice a week). Follow the *easycare* instructions for use.

**Step E6 Enter the *easycare* mode.**

Press and **hold-down** the ON/OFF button (Fig. E6) on the controller until the *easycare* screen appears (Fig. E6b), then release the ON/OFF button. You will hear one "beep" and see the ***easycare* cleaning screen** (Fig. E6c).

Fig. E6


Fig. E6a


Hold ON/OFF button until *easycare* screen appears

Fig. E6b

Release ON/OFF button

Fig. E6c

**Step E7** After backwashing follow the instructions provided in chapter F for disinfecting your nebuliser handset parts and the aerosol head.

## F Disinfecting Your Nebuliser Handset

**Note:** Clean your nebuliser handset before disinfecting your nebuliser handset. The nebuliser handset cannot be disinfected effectively unless it has been cleaned beforehand.

### Method 1 (preferred): Thermal Disinfection (boiling)

Disinfect the nebuliser handset parts by boiling in **distilled water**.

**Step F1** Fill a clean pot with enough **distilled water** to fully cover all parts.


**⚠ CAUTION**

Plastic may melt if it comes into contact with the hot bottom of the pot unless there is enough water in the pot.

**Step F2** Bring the water to a boil.

**Step F3** Add the disassembled parts, including the aerosol head and boil for a **full 5 minutes**.

*Fig. F3*


**Step F4** Remove the parts after the water has cooled and shake the parts to remove excess water.

**Step F5** Air-dry parts on a clean, lint-free towel.

*Fig. F5*

**Step F6** After the parts are completely dry, store them in a clean bag or cover them with a clean, lint-free towel. Reassemble just before treatment.


## Method 2: Steam Disinfection

Steam disinfection can be performed with a commercially available baby bottle disinfecter with an operating time of at least 6 minutes.

**Step F7** Follow the directions for the baby bottle disinfecter and use the full disinfection cycle with the correct volume of water described in the directions.

*Fig. F7*

**Step F8** Remove the parts after the disinfecter has cooled and shake the parts to remove excess water.

**Step F9** Air-dry parts on a clean, lint-free towel.

**Step F10** After the parts are completely dry, store them in a clean bag or cover them with a clean, lint-free towel. Reassemble just before treatment.

### Method 3: Washer-disinfector (Professional Healthcare only)

Use a washer-disinfector according to the requirements set out in the ISO 15883 series of standards.

The effectiveness of this procedure has been validated by using the alkaline-enzymatic cleaning agent neodisher® MediClean forte (Chemische Fabrik Dr. Weigert, Hamburg) as the cleaning agent in conjunction with demineralised water as a neutraliser.

**The nebuliser handset is for single patient use only and must not be shared between patients.** Therefore it is important to group nebuliser handset parts by patient while washing. Position the parts so that they are cleaned thoroughly everywhere.

Use demineralised water and a suitable cleaning agent, following the instructions of the cleaning agent manufacturer. Select the programme for cleaning and subsequent thermal disinfection that achieves at least an A0 value of 3000.

Even if the washer-disinfector is equipped with a drying function, make sure that no residual moisture remains in the nebuliser handset by shaking the parts and allowing them to dry on a clean, lint-free towel. After the parts are completely dry, store them in a clean bag or cover them with a clean, lint-free towel.

Reassemble just before treatment.

## G Cleaning and Disinfecting Your Controller, Connection Cord, Power Supply

**Step G1** If the controller, connection cord or power supply needs cleaning, disconnect the connection cord and power supply from the controller and unplug the power supply from the wall outlet.

**Step G2** **To clean**, wipe down the controller housing, connection cord and power supply with a clean, damp, soft towel.


**⚠ CAUTION**

Do not allow liquid to get inside the controller, because it may damage the electronics.

**Step G3** **To disinfect**, wipe down the controller housing, connection cord and power supply with an alcohol-based disinfecting wipe. To ensure safe application, follow the instructions for use of the disinfectant wipe. Effectiveness has been tested using Bacillol Tissues and Clinell Wipes.

**Step G4** Allow the controller, connection cord and power supply to dry completely before use. Always store the parts in a dry and clean location.

Fig. G3

# H Assembling Your Nebuliser System

## ⚠ WARNING


Your nebuliser handset must be cleaned before first use. Follow the instructions provided in chapter E.

**Note:** Inhaling too cold aerosol can irritate the airways. Take the medication out of the refrigerator before you assemble your nebuliser system.

### Step H1 Wash your hands with soap and water.

Dry them well.

Fig. H1

### Step H2 Insert the aerosol head.


#### ⚠ CAUTION

**Do not** touch the silver center part of the aerosol head. Touching it may damage the membrane that generates the aerosol.

Turn the aerosol head so that the writing e.g. "eFlow<sup>®</sup>rapid" is at the top of the aerosol head. Make sure to insert the accordingly labeled aerosol head into the respective nebuliser handset. Using two fingers, carefully grasp the aerosol head by the two flexible plastic tabs (a) on each side of the aerosol head.


Squeeze the two flexible plastic tabs together while inserting the two metal prongs and the two flexible plastic tabs into the medication reservoir (b). You should feel the flexible plastic tabs "grab" as you insert the aerosol head.

Fig. H2


**Step H3**     **Insert the blue inspiratory valve** into the aerosol chamber (a).  
Make sure the flaps (b) are pointed down and seated properly.

Fig. H3


**Step H4**     **Close the nebuliser handset** by pushing the aerosol chamber (a) together with the medication reservoir (b).  
You should hear a slight "snap" as the tab (c) latches onto the aerosol chamber.  
If you do not hear and feel a slight snap, open the nebuliser handset and check that the blue inspiratory valve is seated properly.

Fig. H4

**Step H5**     **Attach mouthpiece** to the nebuliser handset.  
Make sure the blue expiratory valve (a) is pressed into the slot of the mouthpiece (b).  
Push the mouthpiece straight onto the nebuliser handset with the blue expiratory valve on the top of the mouthpiece.


Fig. H5


**Step H6 Attach the connection cord to the nebuliser handset.**

Line-up the bottom of the connector (a) with the bottom of the nebuliser handset (b).

Lift until the tab on the connection cord (c) snaps into the rear underside of the nebuliser handset.

*Fig. H6***Step H7 Attach the connection cord to the controller.**

Attach the rectangular end of the connection cord (a) to the bottom of the controller.

*Fig. H7*

## I Taking Your Treatment

**Step I1**     **Add the drug to the medication reservoir.** Make sure the nebuliser handset is on a flat, stable surface.


**⚠ CAUTION**

**Do not** attempt to inhale essential oils with an eFlow<sup>®</sup>rapid nebuliser handset, as this can irritate the airways (intense coughing).

**⚠ WARNING**

- **Do not** overfill the nebuliser handset. Always make sure the medication is not above the maximum fill level (see upper scale graduation mark on the nebuliser handset).
- To ensure the correct dose, **never** use another drug for a drug-specific nebuliser handset than the one approved, and **do not** use with an eFlow<sup>®</sup>rapid nebuliser handset.

Fig. I1


**Step I2**     **Attach the medication cap** by aligning the tabs on the medication cap with the tab slots on the medication reservoir.

Turn the medication cap clockwise until it stops.

**Note:** As the medication cap is turned, the inner cap (a) will rise until it matches the height of the outer cap.

**Note:** Before use, note that your device must be properly assembled and the medication reservoir is sealed. Use of an improperly assembled device could decrease the effectiveness of your treatment.


Fig. I2


- Step I3** **To get ready, sit in a relaxed, upright position.**  
Place the mouthpiece on top of your bottom lip and tongue.  
Close your lips around the mouthpiece but **do not** cover the blue expiratory valve (a).  
Hold your nebuliser handset level.

Fig. I3

- Step I4** **Start your treatment.** Press the controller ON/OFF button (a). You will hear one "beep" and the eFlow® Technology logo will appear on the display. The device is now operating.

Fig. I4

- Step I5** You will see a timer running down from 3 to 1 before treatment starts (Fig. I5a).  
Aerosol will begin to flow into the aerosol chamber of the nebuliser handset.  
The controller will display the **treatment screen** (Fig. I5b).

Fig. I5a



Fig. I5b


**Step 16 Breathe normally through the mouthpiece.** The mouthpiece should remain in your mouth during treatment. As you exhale you will see aerosol come out of the blue expiratory valve on the mouthpiece. This is normal.

**Do not** breathe through your nose.

**⚠ CAUTION**

Keep the nebuliser handset level or the controller may switch off before your treatment is complete.


Fig. 16

**Step 17 To pause your treatment.** Press the OK button. The pause mode is activated when the controller display switches to the **pause screen** (Fig. 17b). To continue your treatment, press the OK button again. The treatment will resume when it switches back to the treatment screen.

Fig. 17a


Fig. 17b

**Step 18 At the end of treatment, the controller will shut off.**

You will hear two beeps and the dose complete checkmark will appear briefly on the display (Fig. 18).

Fig. 18



**Step I9** Once the treatment is complete, **take your nebuliser handset apart for cleaning** by reversing the steps used to assemble your nebuliser handset. The aerosol head can be removed by squeezing the two flexible tabs (a) on the sides of the aerosol head.

**Note:** For the eFlow<sup>®</sup> *rapid* nebuliser handset about 1 ml of medication cannot be nebulised and remains unused in the medication reservoir. This is intentional and does not represent a malfunction. For safe disposal of leftover medication, follow the information in the package leaflet accompanying the medication.

**Note:** For drug-specific nebuliser handsets, e.g. Tolero<sup>®</sup> for Vantobra<sup>®</sup>, no significant residual volume remains in the medication reservoir. If you see more than a drop in the medication reservoir, reassemble the nebuliser handset and continue treatment.

Fig. I9

**Step I10** Follow the instructions provided in chapter E for cleaning your nebuliser handset.

## J Controller Settings

You can make settings with regards to the illumination of your display as well as the sound indicators.

**Step J1** **Enter the settings mode.** Press and **hold-down** the ON/OFF button (Fig. J1) until you skip the *easyCare* screen.

When the settings screen appears (Fig. J1b), release the ON/OFF button. You will hear one "beep" and see the **backlight default setting screen** (Fig. J1c).

Fig. J1



Fig. J1a


Hold ON/OFF button until settings screen appears

Fig. J1b


Release ON/OFF button

Fig. J1c


**Step J2**     **Change backlight settings.** In the first step, you can make changes to the backlight. By repeated pressing of the ON/OFF button you can choose between the options "Default" (Fig. J2a), "Backlight always ON" (Fig. J2b) and "Backlight always OFF" (Fig. J2c). In the delivered state, the option "Default" is selected: Backlight switches on for 1 second when button is pressed, or the device status changes.

Fig. J2a


Default

Fig. J2b

Backlight always ON


Fig. J2c

Backlight always OFF

**Step J3**     **Apply preferred setting.** Navigate to your preferred setting. Wait until the controller automatically saves the selection after a few seconds. The controller confirms this with a checkmark (Fig. J3).

Fig. J3

**Step J4** **Change sound settings.** In the second step you can change the intensity of the signal tones. By repeated pressing the ON/OFF button you can choose between the options: "Default" (Fig. J4a), "Reduced signal tones" for one single beep for status feedbacks (Fig. J4b) and "All signal tones OFF" (Fig. J4c). In the delivered state, the option "Default" is selected: All signal tones ON.

Fig. J4a


Default

Fig. J4b

Reduced signal tones

Fig. J4c

All signal tones OFF

**Step J5** **Apply preferred setting.** Navigate to your preferred setting. Wait until the controller automatically saves the selection after a few seconds. The controller confirms this with a checkmark (Fig. J5).





Fig. J5



**Step J6** **The controller will shut off. The preferred settings are now saved.**



## K Troubleshooting

### Controller Signal When Operating Properly





#	Controller Screen Display	Display Symbol/ Controller Sound	Condition	Action Required
1.		Timer running down	<b>Normal:</b> The timer starts before treatment.	<b>No action required.</b> The controller is working properly.
2.		Steady (alternating mist)	<b>Normal:</b> Generating aerosol, no errors detected.	<b>No action required.</b> The nebuliser handset is connected and working properly.
3.		Brief two beep tones	<b>Normal:</b> End of dose.	<b>No action required.</b> The treatment is complete (see Notes, Step I9).
4.		Steady	<b>Normal:</b> Controller has paused.	To resume, press the ON/OFF button as indicated on the screen.

#	Controller Screen Display	Display Symbol/ Controller Sound	Condition	Action Required
5.		Flashing alternating signal tone	<b>Normal:</b> The device was turned off.	If the device was switched off manually and medication remains, restart the controller.
6.		Steady	<b>Normal:</b> The <i>easycare</i> mode was selected.	<b>No action required.</b> The <i>easycare</i> is connected and working properly.

### Controller Signal When Failures Occur and Action is Required

#	Controller Screen Display	Display Symbol/ Controller Sound	Condition	Action Required
7.		Flashing alternating signal tone	<ul style="list-style-type: none"> <li>• Bad or missing connection.</li> </ul>	<ul style="list-style-type: none"> <li>• Ensure, that all parts are connected tightly.</li> <li>• Replace aerosol head and/or the connection cord, if defective.</li> <li>• Restart treatment by pressing ON/OFF button.</li> <li>• If the condition persists, contact your service partner.</li> </ul>
8.		Flashing alternating signal tone	<ul style="list-style-type: none"> <li>• No drug detected.</li> </ul>	<ul style="list-style-type: none"> <li>• Add drug, if necessary.</li> <li>• Hold the nebuliser handset level and restart treatment by pressing the ON/OFF button.</li> </ul>

#	Controller Screen Display	Display Symbol/ Controller Sound	Condition	Action Required
9.		Flashing alternating signal tone	<ul style="list-style-type: none"> <li>The device has reached 20-minute maximum time and will shut off.</li> </ul>	<ul style="list-style-type: none"> <li>Check to see if medication remains in the medication reservoir (see Notes, Step I9).</li> <li>If medication remains, restart the controller.</li> <li>After treatment is complete, clean the aerosol head.</li> </ul>
10.		Flashing alternating signal tone	<ul style="list-style-type: none"> <li>Battery is empty, no nebulisation is possible.</li> </ul>	<ul style="list-style-type: none"> <li>Connect power supply to take / resume treatment and keep power supply plugged afterwards to charge battery.</li> </ul>
11.		Flashing alternating signal tone	<ul style="list-style-type: none"> <li>Battery pack is defective.</li> </ul>	<ul style="list-style-type: none"> <li>Switch to power supply to continue treatment.</li> <li>Contact your service partner.</li> </ul>
12.		Flashing alternating signal tone	<ul style="list-style-type: none"> <li>Battery pack is overheated.</li> </ul>	<ul style="list-style-type: none"> <li>Wait until the device cools down or switch to power supply to continue treatment.</li> </ul>

## Other Failures and Faults Which Require Action

#	Fault/Condition	Possible Cause/Solution
13.	The controller will not power up (no beep).	<ul style="list-style-type: none"> <li>- Connect power supply to charge batteries.</li> <li>- Check to see if the power supply is plugged into the controller and the wall outlet.</li> </ul>
14.	The device switches off prematurely during operation (see screen 3.), even though there is still medication in the medication reservoir.	<ul style="list-style-type: none"> <li>- The nebuliser handset was not being held level. Hold the nebuliser handset level and press the ON/OFF button again to resume treatment.</li> </ul>
15.	The device will not switch off automatically, even though there is no medication left in medication reservoir.	<ul style="list-style-type: none"> <li>- The treatment session has been completed successfully. To manually switch off the device, press the ON/OFF button.</li> <li>- The aerosol head may be dirty, clean the aerosol head.</li> <li>- If the condition persists, contact the your service partner.</li> </ul>
16.	Nebuliser handset leaks.	<ul style="list-style-type: none"> <li>- The handset is not assembled properly. Confirm the aerosol head is placed correctly in the nebuliser handset and that the blue inspiratory valve is placed correctly in the nebuliser handset with both flaps pointing into the nebuliser handset chamber.</li> <li>- It is normal to have some liquid condense in the aerosol chamber and mouthpiece. The amount depends on the volume of the medicine and on your breathing pattern.</li> <li>- If the condition persists, contact your service partner.</li> </ul>
17.	Longer than normal nebulisation time.	<ul style="list-style-type: none"> <li>- The aerosol head was not cleaned immediately after it was used, clean the aerosol head.</li> <li>- In addition, clean the membrane of the aerosol head using the <i>easycare</i> cleaning aid, if available in your country/sales region.</li> <li>- The aerosol head was used longer or more often than planned. For best performance, replace the aerosol head regularly.</li> <li>- If the condition persists, contact your service partner.</li> </ul>
18.	The medication cap cannot be removed from the medication reservoir after the inhalation treatment.	<ul style="list-style-type: none"> <li>- There is too much negative pressure in the medication reservoir after inhalation treatment. Open the aerosol chamber and gently pull the aerosol head so that air is able to flow in. The medication cap can now be removed easily.</li> </ul>

## L Terms and Conditions of Warranty

PARI guarantees to you that if used for its intended purpose your nebuliser system will be free from defects in material or workmanship caused by the manufacturing process for the warranty period of two years.

However this warranty does not cover wearing parts, that is to say device parts that are exposed to normal wear, such as the nebuliser handset and aerosol head.

The warranty shall no longer be valid if:

- the device has been operated or used in a manner inconsistent with the guidelines in the instructions for use,
- damage is evident that is attributable to external factors such as water, fire, lightning strike or similar,
- the damage has been caused by transporting the device incorrectly, or resulted from a fall,
- the device has been handled or serviced incorrectly,
- the serial number on the device has been changed, removed or rendered illegible,
- repairs, adaptations or modifications have been made to the device by persons not authorized by PARI.

If, exceptionally, a defect does occur, PARI will replace the device. The device may be replaced with an identically or at least comparably equipped model.

Replacement of the device shall not serve as grounds for new claims under the warranty. All exchanged used devices or parts will become the property of PARI.

All other claims under the warranty are excluded to the extent permitted by law. In the event of a claim under the warranty, please contact your local service partner.

The original purchase receipt from the dealer must be presented as proof of warranty and ownership. This warranty applies solely for the initial purchaser of the device.

## M Spare Parts and Accessories

Description	Item number
eFlow <sup>®</sup> <i>rapid</i> nebuliser handset (incl. aerosol head)	678G8222
eFlow <sup>®</sup> <i>rapid</i> Aerosol Head	678B2620
Mouthpiece	078G3601
Connection cord	178G6010
Hard shell carrying case	078E8010
PARI Filter valve set <sup>a)</sup>	041G0500
PARI SMARTMASK (adults) <sup>a)</sup>	041G0730
SMARTMASK Kids (children aged 2 and older) <sup>a)</sup>	078G5000
SMARTMASK Baby (babies weighing 2.5 kg and more) <sup>a)</sup>	078G5026
<i>easycare</i> cleaning aid for aerosol heads <sup>a)</sup>	078G6100
Power supply	078B7116

a) The availability of accessories may differ between countries/sales regions.

## N Technical Data

Weight: controller and connection cord	210 g / 8 oz
Dimensions (LxWxH): controller (approx.)	16 x 7 x 2 cm / 6.3 x 2.8 x 0.8 in
Weight: nebuliser handset	55 g / 1.9 oz
Dimensions (LxWxH): nebuliser handset (approx.)	14.5 x 5.0 x 6.3 cm / 5.7 x 2 x 2.5 in
PARI power supply (REF 078B7116 for international power supply)	Input: 100 - 240 V, 50-60 Hz  Output: 5 V  Total power consumption < 3.5 Watt
Sound pressure level (nebuliser handset)	< 33 dB (A)

### eFlow<sup>®</sup>rapid Performance Characteristics

Min. fill volume	2.0 ml						
Max. fill volume	6.0 ml						
MMAD <sup>a)</sup> (Mass Median Aerodynamic Diameter)	4.6 µm						
GSD <sup>a)</sup> (Geometric Standard Deviation)	1.8						
Respirable fraction <sup>a)</sup>	<table border="0"> <tr> <td>&lt; 5 µm</td> <td>55%</td> </tr> <tr> <td>&lt; 2 µm</td> <td>6.5%</td> </tr> <tr> <td>2 - 5 µm</td> <td>48.5%</td> </tr> </table>	< 5 µm	55%	< 2 µm	6.5%	2 - 5 µm	48.5%
< 5 µm	55%						
< 2 µm	6.5%						
2 - 5 µm	48.5%						
Aerosol output <sup>b)</sup>	0.62 g						
Aerosol output rate <sup>b)</sup>	0.5 g/min						

Aerosol data according to ISO 27427; nebulised medium: Salbutamol 2.5 ml. Measured values may vary depending on the medication used.

Recommended fill volume: see package leaflet for the medication

a) Measurement with Next Generation Pharmaceutical Impactor (NGI) at 23 °C and 50% relative humidity. Inspiratory flow: 15 l/min.

b) Measurement with breathing simulator at 23 °C and 50% relative humidity. 500 ml breathing volume, breathing frequency 15 cycles/minute, sinusoidal breathing pattern, inhalation/exhalation ratio 1:1 (with adults, may vary in children).

For performance characteristics for a nebuliser handset specifically approved for and/or supplied with a drug refer to the information about the drug.

**Transport/Storage (unpacked product, between uses)**

Temperature	-25 to 70 °C/-13 to 158 °F
Humidity	0% to 93%
Air Pressure	500 to 1060 hPa/7 to 15 PSI

**Operating Conditions**

Temperature	5 to 40 °C/41 to 104 °F
Relative Humidity (non-condensing)	15% to 93%
Air Pressure	700 to 1060 hPa/10 to 15 PSI

**Batteries**

PARI rechargeable battery-pack (NiMH)	3.6 V (nominal), 2000 mAh
Runtime	About 2 hours of continuous operation

**Device Classification According to IEC 60601-1 & IEC 60601-1-2**

Type of electric shock protection (power supply)	Protection Class II
Degree of protection from electric shock (nebuliser handset)	Type BF
Degree of protection against water ingress per IEC 60529 (controller)	IP 22
Degree of protection when in the presence of flammable mixtures	No protection
Operating mode	Continuous Operation
Environment	Home use & professional healthcare settings

Electromagnetic compatibility data is available in table form upon request from PARI or on the Internet at <https://www.pari.com/fileadmin/Electromagnetic-compatibility-4.pdf>.

## Device materials

This product consists mainly of the following materials:

- Controller housing: Acrylonitrile butadiene styrene copolymer, thermoplastic elastomers, polyethylene
  - Nebuliser handset including aerosol head: Polypropylene, thermoplastic elastomers, polyoxymethylene, silicone rubber, stainless steel
  - Mouthpiece: Polypropylene, silicone rubber
- Not made with natural rubber (latex).

## Lifetime

With appropriate care, the eFlow<sup>®</sup>*rapid* Aerosol Head can be used between 3 and 6 months and the nebuliser handset body for up to 1 year.

However for best performance it may be required to replace the nebuliser handset including aerosol head more often. Follow the instructions of the nebuliser handset. The controller, connection cord and power supply should last 3 years with appropriate care (which is roughly 150 disinfection cycles).

## Disposal

The nebuliser system parts and batteries must be disposed in accordance with local (state, county or municipal) regulations.

## O Symbols

	Medical device
	Single patient - multiple use
	Consult the instructions for use
	Follow the instructions for use
	Manufacturer
	Date of manufacture
	Lot number
	Item number
	Serial number
	Unique Device Identifier
	CE marking: This product complies with the European requirements for medical devices (EU 2017/745).
	Direct current
	Alternating current
	Protection class of the applied part: Type BF
	Minimum and maximum ambient temperature
	Minimum and maximum humidity
	Minimum and maximum air pressure
	Protection class II equipment

**IP22**

The device is protected against solid foreign objects of 12,5 mm Ø and greater and against vertically falling water drops when controller tilted up to 15°.


Do not dispose in domestic waste. This product must be disposed in accordance with the applicable (national, county or municipal) regulations.


PARI reserves the right to make technical changes. Information as of: 2023-09

# Quick Instructions

1. Please read the entire Instructions for Use taking note of all warnings and cautions.
2. Find these parts (for detailed instruction see page 8).


eBase<sup>®</sup>  
Controller


Connection  
Cord

Aerosol  
Head

Nebuliser  
Handset Body

Inspiratory  
Valve


Mouthpiece


Medication  
Cap

## Nebuliser Handset


3. Clean your nebuliser handset parts before you use them for the first time and immediately after each use (for detailed instruction see page 12).
4. Disinfect your nebuliser handset parts after each treatment day (for detailed instruction see page 15).

1  
Rinse each nebuliser handset part for about 5 seconds.



2  
Soak parts for 5 minutes in warm soapy water and then rinse thoroughly.


3  
Disinfect by boiling for 5 minutes in distilled water or use one of the other methods, then air-dry.

# Quick Instructions


## 5. Assemble your nebuliser handset (for detailed instruction see page 19).

Insert the aerosol head and then insert the blue inspiratory valve.



Attach the mouthpiece.



Attach the connection cord to the nebuliser handset and to the controller.


## 6. Take your treatment (for detailed instruction see page 22).

Add the drug to the medication reservoir and attach the medication cap.

Place the nebuliser handset in your mouth and then start the controller.

The controller shows the treatment screen while nebulisation is ongoing.

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