
eFlow[®] rapid nebuliser system with eTrack[®] Controller

en Instructions for use

da Brugsanvisning

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no Bruksanvisning

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sv Bruksanvisning

Patented. See www.pari.com/ip

©2024 PARI Pharma GmbH. All rights reserved.

PARI Connect® is a registered trademark of PARI Medical Holding GmbH. The PARI Pharma logo is a registered trademark of PARI GmbH. eTrack®, eFlow®, Tolero® and Vantobra® are registered trademarks of PARI Pharma GmbH. TouchSpray made under licence from The Technology Partnership plc. The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc.

en	INSTRUCTIONS FOR USE	3
da	BRUGSANVISNING	47
de	GEBRAUCHSANWEISUNG	91
es	INSTRUCCIONES DE USO	135
fi	KÄYTTÖOHJE	179
fr	NOTICE D'UTILISATION	223
it	ISTRUZIONI PER L'USO	269
nl	GEBRUIKSAANWIJZING	313
no	BRUKSANVISNING	359
pt	MANUAL DE INSTRUÇÕES	403
sv	BRUKSANVISNING	447

en INSTRUCTIONS FOR USE

For eFlow[®]*rapid* nebuliser system with eTrack[®] Controller 178G1110

Table of contents

A	Introduction	4
B	Safety Precautions	6
C	Parts of Your eFlow[®]<i>rapid</i> nebuliser system	8
D	Charging Your Controller	10
E	Initial Set-up for Data Transfer	12
F	Manual Data Transfer	14
G	Cleaning Your Nebuliser Handset	15
H	Disinfecting Your Nebuliser Handset	18
I	Cleaning and Disinfecting Your Controller, Connection Cord, Power Supply	21
J	Assembling Your Nebuliser System	22
K	Taking Your Treatment	25
L	Controller Settings	30
M	Troubleshooting	32
N	Terms and Conditions of Warranty	37
O	Spare Parts and Accessories	38
P	Technical Data	39
Q	Symbols	42
	Quick Instructions	44

A Introduction

The eFlow[®]*rapid* nebuliser system is a quiet, lightweight, battery-operated medical device which converts your drug into an aerosol for inhalation.


The drug is placed into the medication reservoir of your nebuliser handset, which feeds to a membrane containing tiny holes. When the device is turned on, the membrane vibrates forcing the drug through the tiny holes to form a fine aerosol mist, which you breathe into your lungs.

When connected to the PARI Connect[®] App the eTrack[®] Controller records and stores nebulisation data such as the start time and date and its duration when a nebulisation takes place. In order to transmit this data to the app or a central server, the controller is equipped with a Bluetooth[®] and with a Wi-Fi module.

The PARI Connect[®] App offers an overview of your therapy. You can enter your individual therapy plan, set reminders and document your therapies. You can enrich the data over time with your personal wellbeing and notes in the app. With the PARI Connect[®] App you can invite a buddy and share your clinical data with your healthcare professionals.

The PARI Connect[®] App can be downloaded for iOS and Android operating systems from the App Store or Google Play Store.

This "**Instructions For Use**" (IFU) contains information and safety precautions for your nebuliser system. This IFU depicts the eFlow[®]*rapid* nebuliser handset and, if not explicitly described, the instructions must also be followed when using a drug-specific nebuliser handset.

The drug-specific nebuliser handset is made available with a specially approved drug and/or via a specific distribution channel.

Take note of any limitations in the package leaflet accompanying the medication.

Before using your eFlow[®]*rapid* nebuliser system, you must read and understand the entire Instructions For Use, and save it for future reference.

If for any reason you do not understand any part of these instructions, contact your service partner before proceeding with your treatment.

Take special note of all safety precautions marked **Warning or Caution**.

Within the **European Union (EU)**, report serious incidents to the manufacturer and the competent authority.

If your health condition for any reason is not improved, or even worsens, seek professional medical advice.

Indications

The eFlow[®]*rapid* nebuliser system is indicated to administer solutions or suspensions for a nebuliser handset which are approved for the treatment of diseases of the lower airways.

Contraindications

No contraindication is known to the manufacturer of this medical device.

Intended Purpose

The eFlow[®]*rapid* nebuliser system is intended to administer solutions or suspensions for a nebuliser handset as an aerosol for inhalation.

Intended use

The eFlow[®]*rapid* nebuliser system must only be used in accordance with its intended purpose.

Intended Patient Group

The eFlow[®]*rapid* nebuliser system is suitable for all age groups but is limited to the patient group defined for the drug. The nebuliser handset is only suitable for persons who are conscious and breathe spontaneously.

Accessories

For babies and infants who cannot use the mouthpiece to inhale, PARI offers appropriate masks. A mouthpiece with filter reduces the release of the exhaled aerosol into the ambient air. To clean the aerosol head by backwashing the membrane, an *easycare* cleaning aid can be used. Accessories are not available in all countries/sales regions. Contact your service partner if you are unsure whether you should use accessories. When using accessories, follow the instructions that came with the accessories.

B Safety Precautions

Please read the following before using this device.

Warning indicates hazards which, without precautionary measures, can result in serious injury or death.

Caution indicates hazards which, without precautionary measures, can result in minor or moderate injury or impair treatment; or cause damage to the device.

WARNING

To reduce the risk of serious injury:

- 1. Do not** use a damaged power supply or handle a power supply with wet hands.
- 2. Do not** use extension cords where they can fall into a bathtub, sink, or other liquid.
- Also, **do not** use the nebuliser system where it can fall into water, e.g. bathtub.
- 4. Do not** reach for the device if it has fallen into water or other liquid. Unplug immediately from the wall outlet. Retrieve the device **only** after it has been unplugged.
- Improper use may cause damage to the product and lead to injuries and illness. Follow these Instructions for Use and all safety instructions contained herein and then save them for future reference.
- To ensure the correct dose, **never** use another drug for a drug-specific nebuliser handset than the one approved, and **do not** use with an eFlow[®]rapid nebuliser handset.
- In order to avoid infections, the nebuliser system is intended for single patient use only. **Do not** share it with others.
- To avoid infections, clean the nebuliser handset including the aerosol head before they are used for the first time.
- To keep the membrane from clogging and to prevent infections, you must clean the nebuliser handset including the aerosol head immediately after each treatment.
- To prevent possible infections, you must also disinfect the nebuliser handset and aerosol head after each treatment day.
- Always unplug the device right after using and before cleaning.
- To avoid the chance of fire, **do not** use where oxygen or combustible gases are being used such as an oxygen tent.
- Use only original PARI parts and accessories. Otherwise this can lead to malfunction including interference of nearby devices.
- The device may be operated in a car, train or airplane. To minimize possible interference issues, please note the following:
 - It must only be used in passenger areas of trains and airplanes.
 - Do not** use near the control systems of airplanes or trains.
 - It must be operated using battery power while in a car.
- 15. Do not** use the device while driving.
- 16. Do not** place the controller next to or on top of other medical devices during operation unless both devices are monitored constantly to make sure both are operating properly.
- 17. Do not** use in areas exposed to elevated magnetic or electrical radiation such as an MRI scanner or high frequency surgical equipment in hospital.
- Close supervision is necessary when this product is used by or near children or the physically or mentally impaired persons.
- Small or broken parts can be swallowed or block the airways and the connection cord poses a strangulation hazard.
- 20. Never** operate the nebuliser system if parts or cables appear damaged, deformed or heavily discolored.

⚠ CAUTION


To minimize risk of damaging the device or impairing treatment:

1. The controller electronics can be damaged if liquid gets into the controller: **Do not** use the device if liquid gets inside the controller.
2. **Do not** disassemble the controller or modify the equipment at any time. There are no serviceable parts.
3. **Do not** allow pets near the nebuliser system.
4. Portable wireless communication devices (such as cell phones or external antennas) may reduce the performance of the nebuliser system. Keep a minimum distance of 30 cm or 12 inches of these devices and the nebuliser system.
5. Anti-theft systems and Radio Frequency Identification (RFID) readers may reduce performance of the nebuliser system: **Do not** use near entrances of stores, libraries, and hospitals.


C Parts of Your eFlow[®]rapid nebuliser system

Make sure you have the parts shown below. Several parts are delivered in a carrying case which can also be used to transport the system. The eFlow[®]rapid nebuliser system contains an *easycare* cleaning aid.


Nebuliser Handset (2x)


Aerosol Head


Nebuliser Handset Body


Medication Cap


Blue Inspiratory Valve


Mouthpiece with blue Expiratory Valve

eTrack[®] Controller


Connection Cord

Power Supply
(EU shown)

General Information About the Nebuliser Handset and Controller

The nebuliser handset body is pre-assembled when delivered and contains the blue inspiratory valve.


You can open the nebuliser handset body by gently pulling up on the tab (c).

The medication reservoir (a) and the aerosol chamber (b) are attached using a "hinge" and must not be detached from one another.

⚠ CAUTION

Do not touch the silver center part of the aerosol head. Touching it may damage the membrane that generates the aerosol.

There are 5 buttons on your eTrack[®] Controller. The ON/OFF button (a) in the top right corner of the device turns the controller on and off. The "left" (◀) and "right" (▶) buttons are used to navigate between different menu items. The BACK button is used to navigate to the previous menu level. The OK button is used to choose the menu item seen on the display.


D Charging Your Controller

The controller runs with an integrated battery pack and can be recharged at any time.

Using the power supply: The power supply will automatically adjust to the incoming voltage and will charge the controller's integrated battery pack. It can be used worldwide, but may require changing the country-specific adapter.

Step D1 To change the country-specific adapter (a), press the clip (b) and slide the country-specific adapter upwards. Slide on the new country specific adapter until the clip connects.


Fig. D1

Step D2 Plug the power supply into the controller.


The power cord socket is located on the top of the controller. Push the end of the power supply cord into the power cord socket (a). Do not attempt to insert the power cord into the bottom of the controller.

Fig. D2

Step D3 Plug the power supply into the wall outlet.

The power supply will charge the batteries in the controller. As soon as the charging process is completed, you will see a fully charged battery.

Fig. D3

Note: To preserve battery life, it is recommended to:

- charge battery only when nearly empty.
- unplug the power supply once the battery is fully charged.
- avoid higher temperatures than room temperature during charging.


E Initial Set-up for Data Transfer

In order for the eTrack[®] Controller to store data and exchange with the PARI Connect[®] App, the controller and the app must be connected to each other via Bluetooth. Then you can connect the controller to your Wi-Fi access point for automatic data transfer.

To connect the controller to your Wi-Fi access point, please have your Wi-Fi password ready in advance. The password needs to be entered into the PARI Connect[®] App.

Step E1 **Turn on the controller.** Press the ON/OFF button (a). You will hear one "beep" and the eFlow[®] Technology logo will appear on the display. The device is now operating and is ready to connect.

Fig. E1

Step E2 **Open your PARI Connect[®] App.** The app guides you through the registration and connection process.

Step E3 Enter PIN when prompted by the app to connect the controller to the PARI Connect® App. The PIN will be shown automatically on the controller display during setup.

Follow the instructions in the app to finalize connection with Bluetooth and proceed with connecting the controller to your Wi-Fi access point.

Alternatively you can connect the controller to your Wi-Fi access point by navigating to the "More" section in the app. Choose "Device Management" and follow the instructions there.

Fig. E3

Note: After successful initialization of your Wi-Fi connection, Bluetooth is deactivated on the controller and nebulisation data is transmitted automatically, if the controller is within the Wi-Fi access point. If data transmission is not possible (e.g. if you are not at home or traveling or due to a connection failure), your data is stored in the controller and sent when communication is possible.

F Manual Data Transfer

If you want to transmit your data stored on the eTrack[®] Controller manually via **Wi-Fi**, do the following:

Step F1 Turn on the controller.

Press the ON/OFF button. You will hear one "beep" and the eFlow[®] Technology logo will appear on the display.

Step F2 Navigate to the manual data transfer screen

(Fig. F2a) and confirm by pressing the OK button. You know the data transfer was successful if you see the communication screen (Fig. F2b) followed by the check screen (Fig. F2c).

Fig. F2a


Fig. F2b



Fig. F2c

If you want to transmit your data manually via **Bluetooth**, turn on the controller and initiate manual data transfer in the app as described there.

G Cleaning Your Nebuliser Handset

⚠ WARNING

Your new nebuliser handset is not ready for use. You must clean your nebuliser handset before you use it for the first time and after every use. You must also disinfect your nebuliser handset at the end of the treatment day.

⚠ CAUTION

Do not touch the silver center part of the aerosol head. Touching it may damage the membrane that generates the aerosol.

Step G1 **Disassemble the nebuliser handset and rinse the aerosol head** for about 5 seconds on each side under warm running tap water (quality of drinking water; about 37 °C or 98 °F, warm to the touch).

Note: Thorough rinsing of both sides of the aerosol head helps to prevent clogging. Clogging may increase nebulisation time. It is important to clean the aerosol head immediately after each use.

Step G2 **Rinse the opened nebuliser handset body (a), the mouthpiece (b), the blue inspiratory valve (c) and the medication cap (d)** under warm running tap water for about 5 seconds to remove medication residue.


Fig. G1



Fig. G2


Step G3 Add a few drops of **clear liquid dish soap** to warm (about 37 °C or 98 °F) tap water in a clean tub or bowl and **soak all parts of the nebuliser handset for 5 minutes.**

⚠ CAUTION

The membrane of the aerosol head can be damaged very easily.

- **Do not** attempt to clean the aerosol head with a brush.
- **Do not** attempt to clean or disinfect the aerosol head in a home dishwasher.
- **Do not** attempt to disinfect the aerosol head in a microwave.

Fig. G3


Step G4 **Rinse all parts thoroughly** under warm running tap water and then inspect them to make sure they are clean and not visibly damaged.

If parts do not appear to be clean, soak them for another 5 minutes as described in Step G3.

Do not use dirty or damaged parts.

Fig. G4


Step G5 **Shake** the parts to remove excess water. **Air-dry parts** on a clean, lint-free towel.

After the parts are completely dry, store them in a clean bag or cover them with a clean, lint-free towel.

Reassemble just before treatment.

Fig. G5


Note: If necessary, clean the membrane of the aerosol head with the *easycare* cleaning aid (not more than twice a week). Follow the *easycare* Instructions for Use.

Step G6 **Navigate to the *easycare* mode** (Fig. G6a) and confirm by pressing the OK button. During backwashing you will see the *easycare* cleaning screen (Fig. G6b).

Fig. G6a

Fig. G6b

Step G7 After backwashing follow the instructions provided in chapter H for disinfecting your nebuliser handset parts and the aerosol head.

H Disinfecting Your Nebuliser Handset

Note: Clean your nebuliser handset before disinfecting your nebuliser handset. The nebuliser handset cannot be disinfected effectively unless it has been cleaned beforehand.

Method 1 (preferred): Thermal Disinfection (boiling)

Disinfect the nebuliser handset parts by boiling in **distilled water**.

Step H1 Fill a clean pot with enough **distilled water** to fully cover all parts.


⚠ CAUTION

Plastic may melt if it comes into contact with the hot bottom of the pot unless there is enough water in the pot.

Step H2 Bring the water to a boil.

Step H3 Add the disassembled parts, including the aerosol head and boil for **a full 5 minutes**.

Fig. H3


Step H4 Remove the parts after the water has cooled and shake the parts to remove excess water.

Step H5 Air-dry parts on a clean, lint-free towel.

Fig. H5

Step H6 After the parts are completely dry, store them in a clean bag or cover them with a clean, lint-free towel. Reassemble just before treatment.


Method 2: Steam Disinfection

Steam disinfection can be performed with a commercially available baby bottle disinfecter with an operating time of at least 6 minutes.

Step H7 Follow the directions for the baby bottle disinfecter and use the full disinfection cycle with the correct volume of water described in the directions.

Fig. H7

Step H8 Remove the parts after the disinfecter has cooled and shake the parts to remove excess water.

Step H9 Air-dry parts on a clean, lint-free towel.

Step H10 After the parts are completely dry, store them in a clean bag or cover them with a clean, lint-free towel. Reassemble just before treatment.

Method 3: Washer-disinfector (Professional Healthcare only)

Use a washer-disinfector according to the requirements set out in the ISO 15883 series of standards.

The effectiveness of this procedure has been validated by using neodisher® MA (Chemische Fabrik Dr. Weigert, Hamburg) as cleaning agent in conjunction with demineralised water as a neutraliser.

The nebuliser handset is for single patient use only and must not be shared between patients. Therefore it is important to group nebuliser handset parts by patient while washing. Position the parts so that they are cleaned thoroughly everywhere. Use demineralised water and a suitable cleaning agent, following the instructions of the cleaning agent manufacturer. Select the programme for cleaning and subsequent thermal disinfection that achieves at least an A0 value of 3000.

Even if the washer-disinfector is equipped with a drying function, make sure that no residual moisture remains in the nebuliser handset by shaking the parts and allowing them to dry on a clean, lint-free towel. After the parts are completely dry, store them in a clean bag or cover them with a clean, lint-free towel.

Reassemble just before treatment.

I Cleaning and Disinfecting Your Controller, Connection Cord, Power Supply

Step I1 If the controller, connection cord or power supply needs cleaning, disconnect the connection cord and power supply from the controller and unplug the power supply from the wall outlet.

Step I2 **To clean**, wipe down the controller housing, connection cord and power supply with a clean, damp, soft towel.


⚠ CAUTION

Do not allow liquid to get inside the controller, because it may damage the electronics.

Step I3 **To disinfect**, wipe down the controller housing, connection cord and power supply with an alcohol-based disinfecting wipe. To ensure safe application, follow the instructions for use of the disinfectant wipe. Effectiveness has been tested using Bacillol Tissues and Clinell Wipes.

Step I4 Allow the controller, connection cord and power supply to dry completely before use. Always store the parts in a dry and clean location.

Fig. I3

J Assembling Your Nebuliser System

⚠ WARNING


Your nebuliser handset must be cleaned before first use. Follow the instructions provided in chapter G.

Note: Inhaling too cold aerosol can irritate the airways. Take the medication out of the refrigerator before you assemble your nebuliser system.

Step J1 Wash your hands with soap and water.

Dry them well.

Fig. J1


Step J2 Insert the aerosol head.

⚠ CAUTION

Do not touch the silver center part of the aerosol head. Touching it may damage the membrane that generates the aerosol.


Fig. J2

Turn the aerosol head so that the writing e.g. "eFlow[®]rapid" is at the top of the aerosol head.


Using two fingers, carefully grasp the aerosol head by the two flexible plastic tabs (a) on each side of the aerosol head. Make sure to insert the accordingly labeled aerosol head into the respective nebuliser handset.

Squeeze the two flexible plastic tabs together while inserting the two metal prongs and the two flexible plastic tabs into the medication reservoir (b). You should feel the flexible plastic tabs "grab" as you insert the aerosol head.

- Step J3** **Insert the blue inspiratory valve** into the aerosol chamber (a).
Make sure the flaps (b) are pointed down and seated properly.

Fig. J3


- Step J4** **Close the nebuliser handset** by pushing the aerosol chamber (a) together with the medication reservoir (b).

You should hear a slight "snap" as the tab (c) latches onto the aerosol chamber.


If you do not hear and feel a slight snap, open the nebuliser handset and check that the blue inspiratory valve is seated properly.

Fig. J4

- Step J5** **Attach mouthpiece** to the nebuliser handset.
Make sure the blue expiratory valve (a) is pressed into the slot of the mouthpiece (b).
Push the mouthpiece straight onto the nebuliser handset with the blue expiratory valve on the top of the mouthpiece.

Fig. J5


Step J6 Attach the connection cord to the nebuliser handset.

Line-up the bottom of the connector (a) with the bottom of the nebuliser handset (b).

Lift until the tab on the connection cord (c) snaps into the rear underside of the nebuliser handset.


Fig. J6

Step J7 Attach the connection cord to the controller.

Attach the rectangular end of the connection cord (a) to the bottom of the controller.

Fig. J7

K Taking Your Treatment

Step K1 **Add the drug to the medication reservoir.** Make sure the nebuliser handset is on a flat, stable surface.


⚠ CAUTION

Do not attempt to inhale essential oils with an eFlow[®]rapid nebuliser handset, as this can irritate the airways (intense coughing).

⚠ WARNING

- **Do not** overfill the nebuliser handset. Always make sure the medication is not above the maximum fill level (see upper scale graduation mark on the nebuliser handset).
- To ensure the correct dose, **never** use another drug for a drug-specific nebuliser handset than the one approved, and **do not** use with an eFlow[®]rapid nebuliser handset.

Fig. K1


Step K2 **Attach the medication cap** by aligning the tabs on the medication cap with the tab slots on the medication reservoir.

Turn the medication cap clockwise until it stops.

Note: As the medication cap is turned, the inner cap (a) will rise until it matches the height of the outer cap.


Note: Before use, note that your device must be properly assembled and the medication reservoir is sealed. Use of an improperly assembled device could decrease the effectiveness of your treatment.

Fig. K2


- Step K3** **To get ready, sit in a relaxed, upright position.**
Place the mouthpiece on top of your bottom lip and tongue.
Close your lips around the mouthpiece but **do not** cover the blue expiratory valve (a).
Hold your nebuliser handset level.

Fig. K3


- Step K4** **Turn on the controller.** Press the controller ON/OFF button (a). You will hear one "beep" and the eFlow[®] Technology logo will appear on the display. The device is now operating.

Fig. K4


- Step K5** **Select the treatment mode** by pressing the OK button.

Fig. K5

If you have entered and synchronized a therapy plan in your app, **choose a drug** by pressing the "left" (◀) or "right" (▶) button and confirm your choice with the OK button (Fig. K5a).


Fig. K5a

You will see a timer running down from 3 to 1 before treatment starts (Fig. K5b).


Fig. K5b

Aerosol will begin to flow into the aerosol chamber of the nebuliser handset.

The controller will display the **treatment screen** followed by the **drug screen** with date (DD.MM.YYYY) and time (24 hr) (Fig. K5c).

Fig. K5c

Step K6


Breathe normally through the mouthpiece. The mouthpiece should remain in your mouth during treatment. As you exhale you will see aerosol come out of the blue expiratory valve on the mouthpiece. This is normal.

Do not breathe through your nose.

⚠ CAUTION

Keep the nebuliser handset level or the controller may switch off before your treatment is complete.

Fig. K6


Step K7 To pause your treatment.

Press the OK button. The pause mode is activated when the controller display switches to the **pause screen** (Fig. K7b). To continue your treatment, press the OK button again. The treatment will resume when it switches back to the treatment screen.

Fig. K7a

Fig. K7b

**Step K8 At the end of treatment, the controller will shut off.**

- You will hear two beeps and the dose complete checkmark will appear briefly on the display (Fig. K8a).
- The communication screen will appear (Fig. K8b).
- The data transfer checkmark will appear briefly on the display (Fig. K8c) if the controller is within range of your Wi-Fi access point and data transfer was successful.
- Then the controller will automatically shut off.

Fig. K8a


Fig. K8b


Fig. K8c



Step K9 Once the treatment is complete, **take your nebuliser handset apart for cleaning** by reversing the steps used to assemble your nebuliser handset. The aerosol head can be removed by squeezing the two flexible tabs (a) on the sides of the aerosol head.

Note: For the eFlow[®] *rapid* nebuliser handset about 1 ml of medication cannot be nebulised and remains unused in the medication reservoir. This is intentional and does not represent a malfunction. For safe disposal of medication residue, follow the information in the package leaflet accompanying the medication.

Note: For drug-specific nebuliser handsets, e.g. Tolero[®] for Vantobra[®], no significant drug residue remains in the medication reservoir. If you see more than a drop in the medication reservoir, reassemble the nebuliser handset and continue treatment.

Fig. K9


Step K10 Follow the instructions provided in chapter G for cleaning your nebuliser handset.

L Controller Settings


Step L1 Turn on the controller. Press the ON/OFF button. You will hear one "beep" and the eFlow® Technology logo will appear on the display.

Fig. L1

Step L2 Navigate to the settings screen and confirm by pressing the OK button.

Fig. L2

- Step L3** By navigating between different settings options, you can adjust the following:
- Change sound settings (Fig. L3a)
 - All signal tones switched on (default)
 - Reduced signal tones
 - All signal tones switched off
 - Change backlight settings (Fig. L3b)
 - Backlight switches on for 6 seconds when button is pressed, or the device status changes (default)
 - Backlight switched on
 - Backlight switched off
 - Switch on/off the airplane mode (Fig. L3c) in order to (de)activate Wi-Fi and Bluetooth (e.g. in situations when usage of Bluetooth or WiFi is not allowed)
 - Reset to factory settings (Fig. L3d)
- Note:** The last menu item displays the device's current date and time (Fig. L3e).

Fig. L3a


Fig. L3b


Fig. L3c


Fig. L3d



Fig. L3e










- Step L4** Settings can be managed by entering the respective item, choosing your desired configuration by pressing the "right" and "left" buttons and then confirm your selection by pressing the OK button.

Note: By resetting the device to factory settings, all stored Wi-Fi and Bluetooth connections are deleted. All stored personal nebulisation data will remain on the device and is not affected by factory reset. Follow the instructions in chapter E for initiating Wi-Fi connection.






M Troubleshooting




Controller Signal When Operating Properly

#	Controller Screen Display	Display Symbol/ Controller Sound	Condition	Action Required
1.		Timer running down	Normal: The timer starts before treatment.	No action required. The controller is working properly.
2.		Steady (alternating mist)	Normal: Generating aerosol, no errors detected.	No action required. The nebuliser handset is connected and working properly.
3.		Brief two beep tones	Normal: End of dose.	No action required. The medication reservoir should be empty (see Notes, Step K9).
4.		Steady	Normal: Controller has paused.	To resume, press the OK button as indicated on the screen.
5.		Flashing alternating signal tone	Normal: The device was turned off.	If the device was switched off manually and medication remains, restart the controller.

#	Controller Screen Display	Display Symbol/ Controller Sound	Condition	Action Required
6.		Steady	Normal: The <i>easycare</i> mode was selected.	No action required. The <i>easycare</i> is connected and working properly.
7.		Steady (flashing arrow)	Normal: Data transfer via Wi-Fi is ongoing	No action required. The controller is working properly.
8.		Steady (flashing arrow)	Normal: Data transfer via Bluetooth is ongoing	No action required. The controller is working properly.
9.		Brief single beep	Normal: Data transfer via Wi-Fi or Bluetooth completed successfully	No action required. The controller is working properly.

Controller Signal When Failures Occur and Action is Required

#	Controller Screen Display	Display Symbol/ Controller Sound	Condition	Action Required
10.		Flashing alternating signal tone	<ul style="list-style-type: none"> Bad or missing connection. 	<ul style="list-style-type: none"> Ensure, that all parts are connected tightly. Replace aerosol head and/or the connection cord, if defective. Restart treatment by pressing OK button. If the condition persists, contact the service partner.
11.		Flashing alternating signal tone	<ul style="list-style-type: none"> No drug detected. 	<ul style="list-style-type: none"> Add drug, if necessary. Hold the nebuliser handset level and restart treatment by pressing the OK button.
12.		Steady	<ul style="list-style-type: none"> The controller is running, but a wrong date and time are shown on the display. 	<ul style="list-style-type: none"> The controller receives the date and time during synchronization with the app. Connect the controller with the app.
13.		Flashing alternating signal tone	<ul style="list-style-type: none"> The device has reached 20-minute maximum time and will shut off. 	<ul style="list-style-type: none"> Check to see if medication remains in the medication reservoir (see Notes, Step K9). If medication remains, restart the controller. After treatment is complete, clean the aerosol head.
14.		Flashing alternating signal tone	<ul style="list-style-type: none"> Battery is empty, no nebulisation is possible. 	<ul style="list-style-type: none"> Connect power supply to take / resume treatment and keep power supply plugged afterwards to charge battery.

#	Controller Screen Display	Display Symbol/ Controller Sound	Condition	Action Required
15.		Flashing alternating signal tone	<ul style="list-style-type: none"> No data has been transmitted for a period of 7 or more days showing the number of stored nebulisations. 	<ul style="list-style-type: none"> Configure Wi-Fi for your eTrack[®] Controller and keep the controller in the range of your Wi-Fi access point for a longer period. The controller will try to connect regularly and transfer data automatically. Perform manual data transfer to your app.
16.		Flashing alternating signal tone	<ul style="list-style-type: none"> Data transfer failed showing the number of stored nebulisations. 	<ul style="list-style-type: none"> Perform your next therapy within the range of your Wi-Fi access point or keep the eTrack[®] Controller near for several hours in order to allow automated data transfer. Perform manual data transfer to your app.
17.		Flashing alternating signal tone	<ul style="list-style-type: none"> Battery pack is defective. 	<ul style="list-style-type: none"> Switch to power supply to continue treatment. Contact your service partner.
18.		Flashing alternating signal tone	<ul style="list-style-type: none"> Battery pack is overheated. 	<ul style="list-style-type: none"> Wait until the device cools down or switch to power supply to continue treatment.

Other Failures and Faults Which Require Action

#	Fault/Condition	Possible Cause/Solution
19.	The controller will not power up (no beep).	<ul style="list-style-type: none"> - Connect power supply to charge batteries. - Check to see if the power supply is plugged into the controller and the wall outlet.
20.	The device switches off prematurely during operation (see screen 3.), even though there is still medication in the medication reservoir.	<ul style="list-style-type: none"> - The nebuliser handset was not being held level. Hold the nebuliser handset level and press the ON/OFF button again to resume treatment.
21.	The device will not switch off automatically, even though there is no medication left in medication reservoir.	<ul style="list-style-type: none"> - The treatment session has been completed successfully. To manually switch off the device, press the ON/OFF button. - The aerosol head may be dirty, clean the aerosol head. - If the condition persists, contact the service partner.
22.	Nebuliser handset leaks.	<ul style="list-style-type: none"> - The handset is not assembled properly. Confirm the aerosol head is placed correctly in the nebuliser handset and that the blue inspiratory valve is placed correctly in the nebuliser handset with both flaps pointing into the nebuliser handset chamber. - It is normal to have some liquid condense in the aerosol chamber and mouthpiece. The amount depends on the fill volume and on your breathing pattern. - If the condition persists, contact the service partner.
23.	Longer than normal nebulisation time.	<ul style="list-style-type: none"> - The aerosol head was not cleaned immediately after it was used, clean the aerosol head. - In addition, clean the membrane of the aerosol head using the <i>easycare</i> cleaning aid, if available in your country/sales region. - The aerosol head was used longer or more often than planned. For best performance, replace the aerosol head regularly. - If the condition persists, contact the service partner.
24.	The medication cap cannot be removed from the medication reservoir after the inhalation treatment.	<ul style="list-style-type: none"> - There is too much negative pressure in the medication reservoir after inhalation treatment. Open the aerosol chamber and gently pull the aerosol head so that air is able to flow in. The medication cap can now be removed easily.

N Terms and Conditions of Warranty

PARI guarantees to you that if used for its intended purpose your nebuliser system will be free from defects in material or workmanship caused by the manufacturing process for the warranty period of two years.

However this warranty does not cover wearing parts, that is to say device parts that are exposed to normal wear, such as the nebuliser handset and aerosol head.

The warranty shall no longer be valid if:

- the device has been operated or used in a manner inconsistent with the guidelines in the instructions for use,
- damage is evident that is attributable to external factors such as water, fire, lightning strike or similar,
- the damage has been caused by transporting the device incorrectly, or resulted from a fall,
- the device has been handled or serviced incorrectly,
- the serial number on the device has been changed, removed or rendered illegible,

- repairs, adaptations or modifications have been made to the device by persons not authorized by PARI. If, exceptionally, a defect does occur, PARI will replace the device. The device may be replaced with an identically or at least comparably equipped model.

Replacement of the device shall not serve as grounds for new claims under the warranty. All exchanged used devices or parts will become the property of PARI.

All other claims under the warranty are excluded to the extent permitted by law. In the event of a claim under the warranty, please contact your local service partner.

The original purchase receipt from the dealer must be presented as proof of warranty and ownership. This warranty applies solely for the initial purchaser of the device.

O Spare Parts and Accessories

Description	Item number
eFlow [®] <i>rapid</i> nebuliser handset (incl. aerosol head)	678G8222
eFlow [®] <i>rapid</i> Aerosol Head	678B2620
Mouthpiece	078G3601
Connection cord	178G6010
Hard shell carrying case	078E8010
PARI Filter valve set ^{a)}	041G0500
PARI SMARTMASK (adults) ^{a)}	041G0730
SMARTMASK Kids (children aged 2 and older) ^{a)}	078G5000
SMARTMASK Baby (babies weighing 2.5 kg and more) ^{a)}	078G5026
<i>easycare</i> cleaning aid for aerosol heads ^{a)}	078G6100
Power supply	078B7116

a) The availability of accessories may differ between countries/sales regions.

P Technical Data

Weight: controller and connection cord	210 g / 8 oz
Dimensions (LxWxH): controller (approx.)	16 x 7 x 2 cm / 6.3 x 2.8 x 0.8 in
Weight: nebuliser handset	55 g / 1.9 oz
Dimensions (LxWxH): nebuliser handset (approx.)	14.5 x 5.0 x 6.3 cm / 5.7 x 2 x 2.5 in
PARI power supply (REF 078B7116 for international power supply)	Input: 100 - 240 V, 50-60 Hz  Output: 5 V  Total power consumption < 3.5 Watt
Sound pressure level (nebuliser handset)	< 33 dB (A)

eFlow[®] rapid Performance Characteristics

Min. fill volume	2.0 ml						
Max. fill volume	6.0 ml						
MMAD ^{a)} (Mass Median Aerodynamic Diameter)	4.6 µm						
GSD ^{a)} (Geometric Standard Deviation)	1.8						
Respirable fraction ^{a)}	<table border="0"> <tr> <td>< 5 µm</td> <td>55%</td> </tr> <tr> <td>< 2 µm</td> <td>6.5%</td> </tr> <tr> <td>2 - 5 µm</td> <td>48.5%</td> </tr> </table>	< 5 µm	55%	< 2 µm	6.5%	2 - 5 µm	48.5%
< 5 µm	55%						
< 2 µm	6.5%						
2 - 5 µm	48.5%						
Aerosol output ^{b)}	0.62 g						
Aerosol output rate ^{b)}	0.5 g/min						

Aerosol data according to ISO 27427; nebulised medium: Salbutamol 2.5 ml. Measured values may vary depending on the medication used.

Recommended fill volume: see package leaflet for the medication

a) Measurement with Next Generation Pharmaceutical Impactor (NGI) at 23 °C and 50% relative humidity. Inspiratory flow: 15 l/min.

b) Measurement with breathing simulator at 23 °C and 50% relative humidity. 500 ml breathing volume, breathing frequency 15 cycles/minute, sinusoidal breathing pattern, inhalation/exhalation ratio 1:1 (with adults, may vary in children).

For performance characteristics for a nebuliser handset specifically approved for and/or supplied with a drug refer to the information about the drug.

Transport/Storage (unpacked product, between uses)

Temperature	-25 to 70 °C/-13 to 158 °F
Humidity	0% to 93%
Air Pressure	500 to 1060 hPa/7 to 15 PSI

Operating Conditions

Temperature	5 to 40 °C/41 to 104 °F
Relative Humidity (non-condensing)	15% to 93%
Air Pressure	700 to 1060 hPa/10 to 15 PSI

Batteries

PARI rechargeable battery-pack (NiMH)	3.6 V (nominal), 2000 mAh
Runtime	About 2 hours of continuous operation

Radio interfaces

Bluetooth	2.4 GHz, typical radiated output power 8 dBm
Wi-Fi	2.4 GHz, typical radiated output power 18 dBm, WPA2
FCC-ID of radio module (Bluetooth/Wi-Fi combination)	XPYNINAW15

Device Classification According to IEC 60601-1 & IEC 60601-1-2

Type of electric shock protection (power supply)	Protection Class II
Degree of protection from electric shock (nebuliser handset)	Type BF
Degree of protection against water ingress per IEC 60529 (controller)	IP 22
Degree of protection when in the presence of flammable mixtures	No protection
Operating mode	Continuous Operation
Environment	Home use and professional healthcare

Electromagnetic compatibility data is available in table form upon request from PARI or on the Internet at <https://www.pari.com/fileadmin/Electromagnetic-compatibility-4.pdf>.

Device materials

This product consists mainly of the following materials:

- Controller housing: Acrylonitrile butadiene styrene, copolymer, thermoplastic elastomers, polyethylene
 - Nebuliser handset including aerosol head: Polypropylene, thermoplastic elastomers, polyoxymethylene, silicone rubber, stainless steel
 - Mouthpiece: Polypropylene, silicone rubber
- Not made with natural rubber (latex).

Lifetime

With appropriate care, the eFlow[®]rapid Aerosol Head can be used between 3 and 6 month and the nebuliser handset body for up to 1 year.

However for best performance it may be required to replace the nebuliser handset including aerosol head more often. Follow the instructions of the nebuliser handset. The controller, connection cord and power supply should last 3 years with appropriate care (which is roughly 150 disinfection cycles).

Disposal

The nebuliser system parts and batteries must be disposed in accordance with local (state, county or municipal) regulations.

Simplified EU Declaration of Conformity (2014/53/EU)

PARI Pharma GmbH hereby declares that the eTrack[®] Controller type of radio equipment complies with Directive 2014/53/EU. The complete text of the EU Declaration of Conformity is available at the following web address:

https://www.pari.com/fileadmin/eTrack_RED_DoC_2.pdf

Q Symbols

	Medical device
	Single patient - multiple use
	Consult the instructions for use
	Follow the instructions for use
	Manufacturer
	Date of manufacture
	Lot number
	Item number
	Serial number
	Unique Device Identifier
	CE marking: This product complies with the European requirements for medical devices (EU 2017/745).
	Direct current
	Alternating current
	Protection class of the applied part: Type BF
	Minimum and maximum ambient temperature
	Minimum and maximum humidity
	Minimum and maximum air pressure
	Protection class II equipment

IP22

The device is protected against solid foreign objects of 12,5 mm Ø and greater and against vertically falling water drops when controller tilted up to 15°.


Do not dispose in domestic waste. This product must be disposed in accordance with the applicable (national, county or municipal) regulations.


PARI reserves the right to make technical changes. Information as of: 2024-05

Quick Instructions

1. Please read the entire Instructions for Use taking note of all warnings and cautions.
2. Find these parts (for detailed instruction see page 8).


eTrack®
Controller


Connection
Cord

Aerosol
Head

Nebuliser
Handset Body

Inspiratory
Valve


Mouthpiece


Medication
Cap

Nebuliser Handset


3. Clean your nebuliser handset parts before you use them for the first time and immediately after each use (for detailed instruction see page 15).
4. Disinfect your nebuliser handset parts after each treatment day (for detailed instruction see page 18).

1
Rinse each nebuliser handset part for about 5 seconds.



2
Soak parts for 5 minutes in warm soapy water and then rinse thoroughly.


3
Disinfect by boiling for 5 minutes in distilled water or use one of the other methods, then air-dry.

Quick Instructions


5. Assemble your nebuliser handset (for detailed instruction see page 22).

Insert the aerosol head and then insert the blue inspiratory valve.



Attach the mouthpiece.



Attach the connection cord to the nebuliser handset and to the controller.

6. Take your treatment (for detailed instruction see page 25).



Add the drug to the medication reservoir and attach the medication cap.

Place the nebuliser handset in your mouth and then start the controller.

The controller shows the treatment screen followed by the drug screen with date (DD.MM.YYYY) and time (24 hr) while nebulisation is ongoing.

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